



... Please fill in completely and keep a copy ...

MediMatch dental laboratory



www.medimatch.co.uk
T: 020 3875 8530

MHRA: 8879 • Damas • DLA member • Unit 2, The Works, Colville Road, London, W3 8BL • United Kingdom

1 Dr. _____ (name) Acc./PO no.: _____
Post code: _____ - _____ Mobile: _____

IMPLANT PRIVATE

2 *No patient personal data*
Your case reference: _____ Age: _____ Initials: ____ / ____

3 **Stage** Always give us **12 full** working days **Delivery Date**

..... U / L Date: ____ / ____ / ____
 U / L Date: ____ / ____ / ____
 U / L Date: ____ / ____ / ____
 U / L Date: ____ / ____ / ____
 U / L Date: ____ / ____ / ____

4 **Shade** - use VITA guide

5 **Indicate implant system + Platform** (please circle)

LAB USE ONLY
for lab use only, real app.

..... - @
 - @

6 System _____ acc. no. _____
Platform _____ on no. _____

For lab use only, Opened by:

1	2
4	3

7 **Restoration**

will be Cement retained
 Screw retained
 Crown
 Bridge
 Full-ceramic
 zirconia build up (aesthetic ++)
 zirconia full contour (aesthetic --)
 Porcelain Bonded
 non precious (Co-Cr)
 semi precious (Pd)
 precious, gold (Au)
 Temporary Crown
 Acrylic
 Composite Bonded

8 **Abutment:**


CAD/CAM abutment made **MediMatch*** * please check availability and pricing
 ordered from supplier

9 **I have sent:** (always send bite on implant/preparation)

	enclosed	please order	amount of		enclosed	please order
Abutment	<input type="checkbox"/>	<input type="checkbox"/>	_____	Final screw:		
Lab screw	<input type="checkbox"/>	<input type="checkbox"/>	_____	Ti	<input type="checkbox"/>	<input type="checkbox"/>
Analogue	<input type="checkbox"/>	<input type="checkbox"/>	_____	Gold	<input type="checkbox"/>	<input type="checkbox"/>

I have disinfected the impression with: _____
by: _____

Implants only



Please always check your impressions, poor impressions may be liable to a minimum handling charge of £10
 This case is a remake case
 I have enclosed new or old components _____

Custom Made Device | Supplied in an unsterilized state | Terms and conditions apply | Visit medimatch.co.uk


Dental appliance information and delivery note: Dental product designed to satisfy the information, properties and detail of what has been prescribed by the above dentist. Purely for use for the patient described with the above reference. The product meets the general safety and performance requirements entailed in the Annex I and Annex VIII of the Medical Devices Directive and the United Kingdom Medical Devices Regulations act SI 2002 No. 618. This product may have been produced in either or both UK and China by MediMatch.

Instructions for use, storing & handling: It is highly recommended that the product is stored in a clean and safe environment if not used immediately. It is also advised that there must be no contact with materials, liquids or acids that could cause disfiguration or damage to the product. The product should not be subjected to extreme heat. Where applicable, you should take care not to damage the dental piece(s) when removing from its model. **Prescriber Feedback:** To enable our laboratory to comply with the medical devices regulations for post market surveillance, please inform us of any feedback or issues regarding the enclosed device as soon as possible.

THIS DENTAL APPLIANCE IS SUPPLIED IN AN UNSTERILIZED STATE

IMPLANT PRIVATE

STANDARD ABUTMENTS
(Fixed shape and angulation)



CUSTOM ABUTMENTS

How to choose?

Standard abutment Pros	Standard abutment Cons	Custom abutment Pros	Custom abutment Cons
May reduce the cost of the restoration.	Margin placement of the crown is predetermined and cannot be precisely controlled to achieve the best results possible.	It will fit perfectly on the implant, will support the gum tissue, and is milled to your precise specifications.	May be higher in cost.
Standard abutments have a standard range, and each can be adjusted. There are different types of stock abutments that will work with most cases.	The angled abutments never seem to quite match the actual angulation required to yield a proper profile, and require lots of adjustment.	May be used with both bone-level or tissue-level implants.	We cannot mill all implant systems, there is a limited list of implant systems that we can offer as cad cam.
Original implant connection.	It is not ideal for the aesthetic zone, it may be difficult for the lab technician to achieve an optimal emergence profile and have the gingival tissue be supported properly with respect to the adjacent teeth.	Allows for an ideal crown preparation with proper angulation corrections and a more lifelike emergence profile. The margin is able to rest at the most aesthetic level with the tissue, making it much easier for the patient to cleanse and maintain healthy gingiva around the implant.	Implant connection is a replica/compatible connection and is not the original connection.

Custom Made Device | Supplied in an unsterilized state | Terms and conditions apply | Visit medimatch.co.uk

Dental appliance information and delivery note: Dental product designed to satisfy the information, properties and detail of what has been prescribed by the above dentist. Purely for use for the patient described with the above reference. The product meets the general safety and performance requirements entailed in the Annex I and Annex VIII of the Medical Devices Directive and the United Kingdom Medical Devices Regulations act SI 2002 No. 618. This product may have been produced in either or both UK and China by MediMatch.
Instructions for use, storing & handling: It is highly recommended that the product is stored in a clean and safe environment if not used immediately. It is also advised that there must be no contact with materials, liquids or acids that could cause disfiguration or damage to the product. The product should not be subjected to extreme heat. Where applicable, you should take care not to damage the dental piece(s) when removing from its model. **Prescriber Feedback:** To enable our laboratory to comply with the medical devices regulations for post market surveillance, please inform us of any feedback or issues regarding the enclosed device as soon as possible. THIS DENTAL APPLIANCE IS SUPPLIED IN AN UNSTERILIZED STATE